

DSC Cleaner

Dampening System Cleaner Concentrate



Product Usage

Cleaner Concentrate for Fountain Solution Circulations

Features and Benefits

- Intensive cleaning action
dissolves also dried and sludge residues of paper, gum and ink
- Very economical: Concentrate dilution 1:5 up to 1:9 with water depending on the quantity of residues to remove
- Strong antibacteria and antifungal action
- Easy to wash out with water

Packaging

10x1 liter per carton, 10 liter cans

Directions

Dilute Dampening System Cleaner DSC with water ca. 1:5 to 1:9 (depending on the grade of contamination in the system). Fill the system with this solution and pump a couple of hours through the fountain circulation system. Then drain and add water and fountain concentrate again.

Technical Specifications

- form: liquid
- Density: 1,025 g/cm³
- Color: yellow - transparent
- Odour: characteristic
- Water miscible: yes
- PH value: 12
- labelling: Xn Harmful

Storage and Shelf Life

Storage: not under 6°C, not above 40°C
Shelf Life: 12 months
Keep container tightly closed
Keep away from direct sun light