

Rollo Vital Paste/GEL/liquid

Ink Roller Rejuvenator



Product Usage

Ink Roller Cleaner and Rejuvenator

Features and Benefits

- Dissolves and removes dried ink, glaze and water soluble residues (paper, dust, etc.) thoroughly and quick
- cleans deeply but gentle
- with highly efficient solvents
- contains rubber rejuvenators that give back to the roller its former smooth surface

Packaging

ROLLO-VITAL Liquid: 10 x 1 l Bottles
ROLLO-Vital Paste: 8 x 700 ml Cans
ROLLO-Vital Gel: 12 x 500 ml Tubes

Directions

Wash down inking unit, put blade in off position, apply Rollo-Vital onto the rollers and let press run for ca. 5 minutes until residues are dissolved, wash down again with a strong wash. Use 1-2 times in the week or for quick colour change dark to light

Technical Specifications

- Density: 1,1-1,3
- Flash Point: -
- Color: yellowish-brown
- Odour: mild soapy
- Water miscible: NO
- Evaporation: NO

Storage and Shelf Life

Store in a cool and dry place

Expire time: 12 months